

# VORTRAN<sup>®</sup> MANOMETER

*"Provides Air Pressure Indication"*

VORTRAN Manometer is a disposable, analog manometer with a connecting port option of either a 22mm T-piece or a 15mm vertical port.

- It is MR conditional and can be used in an MRI/CT environment.
- Attaches directly to an automatic resuscitator, such as the VORTRAN GO<sub>2</sub>VENT<sup>®</sup>, or a hyperinflation resuscitation bag.



ORDER No.  
6030



ORDER No.  
6031

## Specifications:

**Range:** 0 to 60 cm-H<sub>2</sub>O

**Connection:** 22mm Tee (Model 6030)

15mm Vertical (Model 6031)

**Accuracy:** ± 2 cm H<sub>2</sub>O from 0 - 20 cm H<sub>2</sub>O

± 3 cm H<sub>2</sub>O from 20 - 40 cm H<sub>2</sub>O

± 5 cm H<sub>2</sub>O from 40 - 60 cm H<sub>2</sub>O